

GB Non-contact thermometer	03
PL Termometr bezdotykowy	08
BG Безконтактен термометър	13
CZ Bezkontaktní teplomer	18
DE Berührungsloses Thermometer	23
EE Kontaktivaba termomeeter	28
ES Termómetro sin contacto	33
GR Ανέπαφο θερμομέτρο	38
HR Beskontaktni termometar	44
HU Érintés nélküli hőmérő	49
LT Nekontaktinis infraraudonųjų termometras	54
LV Bezkontakta termometrs	59
RO Termometru fara contact	64
RU Бесконтактный термометр	69
SK Bezkontaktný teplomer	75
UA Безконтактний термометр	80

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Non-contact infrared thermometer

THE MANUFACTURER RESERVES A RIGHT TO MODIFY THE PRODUCT SPECIFICATION WITHOUT NOTICE.

I. SAFETY PRECAUTIONS

- Follow maintenance recommendations stipulated in this manual.
- The device is intended for purposes specified in this manual only.
- The device may operate at the ambient temperature from 10 °C to 40 °C.
- Keep the device in a dry and clean place.
- Do not expose the thermometer to the impact of current.
- It is forbidden to expose the thermometer to extreme temperatures, i.e. $> 55\text{ °C} < -20\text{ °C}$.
- Do not operate the device if the relative humidity is $> 85\%$.
- The most sensitive part of the product is a protective glass that secures a lens.
- Do not touch the protective glass that secures the lens with your fingers.
- Clean the glass with a cotton swab moistened with a 95-% alcohol.
- Do not expose the thermometer to sunlight or water.
- Do not drop the device.
- Do not use the product if its probe is damaged.
- In case of any problems with the device, contact the dealer.
- Do not repair the thermometer on your own.
- Information on disposal of waste, worn-out devices and accessories has been showed in the user's manual.

II. INTENDED USE

This is an infrared thermometer used to take the temperature from the forehead and ear of children and adults, with no need to touch the body.

III. INTRODUCTION

The non-contact infrared thermometer has been designed on the basis of the latest infrared technology. Since the product assures precise, fast and contactless measurement, it is a perfect tool for safe temperature measurement from the forehead and ear canal. Similar to other thermometers, AET-R842 must be used properly to obtain reliable and stable results. For this reason it is advisable to read this user's manual and safety principles before you start using the thermometer.

IV. SAFETY MEASURES BEFORE USE

- The thermometer has been configured in the production plant. You are not required to set it up before you switch it on.
- To make sure low or high ambient temperatures do not affect its precision, put the product in the place with a room temperature of $15\text{ °C} - 40\text{ °C}$ and keep it there for at least 30 minutes before you start using it to avoid erroneous measurements.
- The person whose temperature is taken must stay at the room temperature for at least 20 minutes to regulate the body temperature.

V. PRINCIPLE OF OPERATION

Regular temperatures as per the measuring method
MEASURING METHOD REGULAR TEMPERATURE °C

EAR 35,8 °C ~ 38 °C

FOREHEAD 35,8 °C ~ 37,8 °C

The human body temperature fluctuates throughout the day. It can be affected by external factors, such as age, sex, type and thickness of skin, etc.

VI. HANDY NOTES ON TEMPERATURE MEASUREMENT

- To obtain precise temperature values, every operator must be properly informed about how to take the temperature using the device and must undergo a practical training.
- Remember not to underestimate the temperature measurement procedures although they seem simple.
- Take the temperature in a neutral context. Before you take the temperature, make sure the patient does not take any dynamic physical activity and the room temperature is moderate.
- While evaluating temperature results, remember about physiological temperature fluctuations: body temperature rises by 0.5 °C between 6:00am and 3:00am. Women have higher body temperature by on average 0.2 °C . Women's body temperature changes depending on their ovulation stage. It grows by 0.5 °C in the second half of the cycle and

early pregnancy phases.

- After intensive physical effort, take a rest for at least 30 minutes before you take the temperature.
- Before you start, wipe the skin dry and draw hair aside.
- Remember to have a 5-second break between subsequent measurements.
- Before use, make sure the probe is clean.

DISPLAY INFORMATION

1. Measuring spot
2. Temperature unit
3. Temperature data
4. Memory symbol
5. Battery symbol

VII. HOW TO USE


- Install batteries.
- Before the first use or after putting batteries in, wait for 10-15 minutes to let the device adapt to the ambient temperature.
- Take the batteries out if the device is not used for a long time.

VIII. HOW TO TAKE THE TEMPERATURE


Press the „SCAN“ key to turn the thermometer on. The screen will activate and the last temperature result will be displayed. When you can feel vibration and see ---°C in the screen, it means that the thermometer is ready for work. There are two modes you can choose: ear and forehead.

• Temperature from ear


To take the temperature from ear, remove the forehead measurement cap, you will see  symbol and a light will go on. Put the probe into the ear canal and press the „SCAN“ key. If the measurement is correct, you will feel the vibration and see temperature result in the screen.

A - Ear temperature – babies below the age of 1: Pull the ear back.


B - Ear temperature – children above the age of 1: Pull the ear and tilt it backwards.

• Forehead temperature

Put the cap onto the thermometer and you will see  symbol. Put the thermometer head against the patient's forehead and press the „SCAN“ key. If the measurement is correct, you will feel vibration and see the temperature result.

C - Point the thermometer at the forehead at a distance of 0-3 cm.

IX. MEASUREMENT TIPS

- If the ambient temperature goes beyond 15–40 °C while you are taking the temperature, the temperature result may not be displayed.
- If the temperature below 32°C is displayed, you will see „Lo” and the thermometer will emit 2 vibrating alerts.
- If the temperature is 37,8°C, you will feel 6 subsequent vibrating alerts.
- If the temperature is above 42,2°C, you will see „Hi” and the thermometer will emit 2 vibrating alerts.

X. Switch between Celsius and Fahrenheit

When the product is off, press and hold the „SCAN” key for about 8 seconds. You will see „---” symbol and a current temperature unit will be displayed. Press the „SCAN” key to select °C or F.

XI. HOW TO CALL THE TEMPERATURE VALUE IN MEMORY MODE

When the product is off, press and hold the „SCAN” key for about 4 seconds to make the thermometer enter the memory mode and display the latest temperature value. Press the „SCAN” key to display another value, from the last to the first. The thermometer has a memory that allows calling values.

XII. BATTERY REPLACEMENT

When a battery symbol is displayed, it means that you need to replace a battery. How to replace batteries: open a battery compartment and replace them with new equivalents, pay attention to proper positioning. Wrong placement of the batteries may lead to damage to the thermometer and invalidation of the guarantee. It is forbidden to use storage batteries. Only non-rechargeable batteries are allowed.

XIII. TECHNICAL SPECIFICATION

Product name	Non-contact thermometer Model: AET-R842
Normal working conditions	Ambient temperature: 10°C ~ 40°C (50°F ~ 104°F) Relative humidity: ≤ 85% Pressure: 700 hPa to 1060 hPa
Storage and transportation conditions	Ambient temperature: -20°C ~ 55°C (-4°F ~ 131°F) Relative humidity: ≤ 95%
Batteries	DC 3V (2 AAA batteries)
Measuring range	32.0°C ~ 34.9°C (89.6°F ~ 94.8°F) ± 0.3°C (±0.6°F) 35.0°C ~ 42.0°C (95.0°F ~ 107.6°F) ± 0.2°C (±0.4°F) 42.1°C ~ 43.0°C (107.8°F ~ 109.4.°F) ± 0.3°C (±0.6°F)
Measuring accuracy (forehead)	≤ 3cm (1.2cala)
Measuring spot	Forehead, ear canal
Unit of measure	Celsius and Fahrenheit
Measuring time	1 s
Automatic shutdown	60 ≤ ± 10 s
Memory	32 results
Dimensions	3,6 x 4,9 x 15,7 cm
Weight	66 g (with no batteries)

XIV. CONTENT OF THE SET

Infrared thermometer, bag, 2 AAA batteries, operating manual.

XV. DESCRIPTION OF SYMBOLS

SYMBOL	REFERENCE
	IEC 60417-5333, Type BF applied part
	IEC 60417-5031 Direct current
	Refer to instruction manual / booklet

SYMBOL	REFERENCE
	DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	Specifies serial number

XV. NOTES

- It is not intended for use in the oxygen-rich environment.
- It is not intended for use with inflammable substances.
- It is not intended for use with inflammable anesthetic products.

XVI. ERROR MESSAGES

MESSAGE	PROBLEM	SOLUTION
Hi	Temperature is higher than 42.2°C (108°F)	Use the thermometer only at the specific temperature. In case the message recurs, contact the seller or customer service department.
Lo	Temperature is lower than 32°C (89.6°F)	Use the thermometer only at the specific temperature. In case the message recurs, contact the seller or customer service department.
Err	System failure	Contact the seller or customer service department.
	Low battery	Replace batteries.
ErH	Ambient measuring temperature is too high	Lower the ambient temperature and keep it at 15°C- 40°C.
ErL	Ambient measuring temperature is too low	Raise the ambient temperature and keep it at 15°C- 40°C.

XVII. EMC DECLARATION

Note: Do not use the device that neighbors other devices or is positioned on them as it may lead to malfunctioning of the product. However, if it is necessary, such devices must be supervised while they are operating to make sure their performance is correct.

Note: Using accessories, converters and cables other than specified or delivered by the device manufacturer may lead to increased electromagnetic wave emission or lower electromagnetic immunity of the device, which eventually may cause it to malfunction.

Note: Do not use mobile devices for radio-wave communication (including peripherals, such as antenna cables or outer antennas) at a distance that is lower than 30 cm from any thermometer element; it applies to cables specified by the manufacturer too. Otherwise the performance of the device may be affected.

All information is used to maintain BASIC SAFETY MEASURES and obtain OPTIMAL EFFICIENCY with regard to electromagnetic interference at the designed period of operation. Mobile radio communication devices may affect operation of the thermometer. While using the product, avoid areas with strong electromagnetic interference, e.g. coming from mobile phones, microwave ovens, etc.

Guidelines and manufacturer's declarations – electromagnetic emission	
Emission test	Compliance
RF CISPR11 emission	Group 1
RF CISPR11 emission	B class
IEC 61000-3-2 harmonic emissions	Not applicable
Voltage fluctuations/flicker emission IEC 61000-3-3	Not applicable

Guidelines and manufacturer's declarations – electromagnetic compatibility		
Electromagnetic compatibility	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
IEC 61000-4-4 fast transient/burst testing	Power lines ± 2 kV Input/output lines ± 1 kV	Not applicable
IEC 61000-4-5 surge testing	Not applicable	Not applicable
IEC 61000-4-11 voltage dips, short interruptions and voltage variations testing	0% 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle I 70% 25/30 cycles One-phase at 0 0% 300 cycles	Not applicable
IEC 61000-4-8 power frequency magnetic field testing	30 A/m 50 Hz/ 60 Hz	30 A/m 50 Hz/ 60 Hz
IEC 61000-4-6 conducted disturbances, induced by radio-frequency field testing	150 kHz to 80 MHz 3 Vrms 6 Vrms (in ISM bands) 80% Am at 1 kHz	Not applicable
IEC 61000-4-3 Radiated electromagnetic field testing at radio frequencies	10 V/m 80 MHz – 2,7 Ghz 80% AM at 1 kHz	10V/m 80 MHz – 2,7 Ghz 80% AM at 1 kHz
Please pay attention to the fact that UT is a voltage in the power supply network before application of the test level.		

Radiated electromagnetic field at radio frequencies IEC 61000-4-3 (research specification for EMC resilience of housing sockets for wireless communication devices utilizing radio frequency waves)	Test frequency (MHz)	Band (MHz)	Service	Modulation	Modulation	Distance (m)	Immunity test levels (V/m)
	385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
	450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sinus.	2	0.3	28
	710 745 780	704-787	Band LTE 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
	810 870 930	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, Band LTE 5	Pulse modulation 18 Hz	2	0.3	28
	1720 1845 1970	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, band LTE 1, 3, 4, 25, UMTS	Pulse modulation 217 Hz	2	0.3	28
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, band 7 LTE	Pulse modulation 217 Hz	2	0.3	29
	5240 5500 5785	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9

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Radiated electromagnetic field at radio frequencies IEC 61000-4-3 (research specification for EMC resilience of housing sockets for wireless communication devices utilizing radio frequency waves)	Test frequency (MHz)	Band (MHz)	Service	Modulation	Modulation	Distance (m)	Immunity test levels (V/m)
	385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
	450	430-470	GMRS 460, FRS 460	FM \pm 5 kHz deviation 1 kHz sinus.	2	0.3	28
	710 745 780	704-787	Band LTE 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
	810 870 930	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, Band LTE 5	Pulse modulation 18 Hz	2	0.3	28
	1720 1845 1970	1700- 1990	GSM 1800, CDMA 1900, GSM 1900, DECT, band LTE 1, 3, 4, 25, UMTS	Pulse modulation 217 Hz	2	0.3	28
	2450	2400- 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, band 7 LTE	Pulse modulation 217 Hz	2	0.3	29
	5240 5500 5785	5100- 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9

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